

Midmark M11 Sterilizer Planned Maintenance Check List

<input type="checkbox"/>	Date:	PM Performed By: _____ Max Temperature Readings: Gauge _____ Actual _____
<input type="checkbox"/>	Replace Door Gasket MIG028-799	Remove the old door gasket. Clean the gasket groove in the door, and chamber face. For easy installation of the new gasket, use the “South-East-West-North” method (refer to instructions included with the PM Kit MIK080-808)
<input type="checkbox"/>	Replace Dam Gasket MIG036-800	Use CAUTION when removing or installing the dam gasket – the door studs can be sharp. To install the door dam, align and place the dam gasket on and over the door studs. The dam must lie flat against the door (refer to the Installation Instructions included with the MIG036-800 Dam Gasket).
<input type="checkbox"/>	Replace Door Springs MIS076-816	Remove the old spring by removing the screw holding it in place. Reuse the screw to install the two (2) new door springs – one on top of the other. This will make the door opening easier.
<input type="checkbox"/>	Replace Filters MIF038-796 or MIF062-798	It is necessary to remove all trays and the rack from the chamber. Clean the chamber then remove the old chamber filter and install the new filter. Depending on when your autoclave was manufactured, you will need to replace either the sintered bronze style Chamber Filter (MIF038-796) or both the mesh style Chamber Filter and Vent Filter (MIF062-798). Chamber Filter (sintered bronze style filter) – MIF038-796: If the machine has only a Chamber Filter Chamber & Vent Filters (mesh style filter) – MIF062-798. If the machine has both a Chamber and a Vent Filter, the Vent Filter serves is located in the vent line to the tee fitting. Reinstall rack and trays.
<input type="checkbox"/>	Sterilizer	After replacing the PM parts: <ul style="list-style-type: none"> <input type="checkbox"/> Check The Reservoir & Clean If Dirt Or Debris Are Present (Note: Debris is the number one reason for Steam Trap/Bellows Failure <input type="checkbox"/> Refill Reservoir With Distilled Water Only <input type="checkbox"/> Conduct A Temperature check Using A Maximum Registering Thermometer (Part #RPT113-838) <input type="checkbox"/> Record The Gauge Temperature & Reading From The Maximum Registering Thermometer <input type="checkbox"/> Conduct A Pressure Check To Identify Any Leaks And Correct Them

Visit www.AllClaveParts.com For Replacement Parts And Assistance
Retain this completed check list for future scheduling and as a record for inspections